
# INSTRUCTIONS FOR USE

Ambu® aScope™ 5 Cysto HD

Ambu

**QUICK GUIDE**

# CONTENTS

# PAGE

English (Instructions For Use).....	4-11
Česky (Návod k použití).....	12-20
Dansk (Brugsanvisning).....	21-29
Deutsch (Bedienungsanleitung).....	30-38
Ελληνικά (δηγίες Χρήσεως).....	39-47
Español (Manual de instrucciones).....	48-56
Suomi (Käyttöohje).....	57-64
Français (Mode d'emploi).....	65-73
Magyar (Használati útmutató).....	74-82
Italiano (Manuale d'uso).....	83-91
日本語 (使用法).....	92-99
Nederlands (Gebruiksaanwijzing).....	100-108
Norsk (Brukerveiledning).....	109-116
Polski (Instrukcja obsługi).....	117-125
Português (Manual de instruções).....	126-134
Slovenčina (Návod na použitie).....	135-143
Svenska (Instruktionshandbok).....	144-152
Türkçe (Kullanım talimatları).....	153-160
中文 (使用说明).....	161-167

## 1. Important information – Read before use

In this document the term aScope 5 Cysto HD refers to Ambu® aScope™ 5 Cysto HD.

Read these Instructions For Use carefully before using aScope 5 Cysto HD. The Instructions For Use may be updated without notice. Copies of the current version are available upon request. Note that these instructions do not explain or discuss clinical procedures.

Before initial use of the endoscope, it is essential for operators to have received sufficient training in clinical endoscopic techniques and to be familiar with the intended use, indications, warnings, cautions, and contraindications mentioned in these instructions.

aScope 5 Cysto HD is not covered by a warranty.

In this document “aScope 5 Cysto HD” refers to instructions that apply to the cysto-nephroscope only and “aScope 5 Cysto HD system” refers to information relevant to aScope 5 Cysto HD, compatible Ambu displaying units and accessories. Unless otherwise specified, endoscope refers to all aScope 5 Cysto HD variants.

### 1.1. Intended use

aScope 5 Cysto HD is a sterile, single-use, flexible cysto-nephroscope intended to be used for endoscopic access to and examination of the lower urinary tract and kidney.

The cysto-nephroscope is intended to provide visualization via a compatible Ambu displaying unit and can be used with endoscopic accessories and instruments.

#### 1.1.1. Intended patient population

aScope 5 Cysto HD is designed for use in adults requiring cystoscopy and/or nephroscopy.

#### 1.1.2. Intended use environment

aScope 5 Cysto HD is designed for use in a hospital environment or medical office environment. Hospital environment includes:

- Out-patient clinics
- Patient wards
- Operating Rooms (OR)
- Intensive Care Units (ICU)
- Emergency Rooms (ER)

### 1.2. Indications for use

aScope 5 Cysto HD is used in patients with disease in the lower urinary tract or kidney who require visualization and examination with a flexible cysto-nephroscope and the use of endoscopic accessories and instruments.

### 1.3. Intended user

aScope 5 Cysto HD is for use by healthcare professionals: Medical doctors, urologists, surgeons, trained within cystoscopy procedure and/or nephroscopy procedure, as well as nurses under medical responsibility trained within cystoscopy.

### 1.4. Contraindications

No contraindications are identified for aScope 5 Cysto HD.

### 1.5. Clinical benefits

Together with the compatible Ambu displaying unit, aScope 5 Cysto HD provides endoscopic access and visualization, enabling examination and procedures in the lower urinary tract and kidney.

### 1.6. Warnings and cautions

#### WARNINGS

1. Only to be used by healthcare professionals trained in clinical endoscopic techniques and procedures specific to urinary tract endoscopy and in accordance with the intended use of the device. Failure to comply with this may cause patient injury.

2. The user must exercise professional judgement when deciding whether a cystoscopy or nephroscopy procedure will be beneficial and necessary for high risk patients.
3. Do not use the aScope 5 Cysto HD if the inspection and preparation, of the product fails as it can cause patient injury.
4. Do not attempt to reuse, reprocess or sterilize the aScope 5 Cysto HD as it is a single use device. Reuse, reprocessing or resterilization of the product may cause structural and functional damage which may cause patient injury including but not limited to mucosal abrasion or bleeding. Reuse, reprocessing or resterilization of the product may also cause contamination leading to infections or sepsis. Cleaning residues left on the product may cause allergic reactions.
5. The distal end of the aScope 5 Cysto HD may get warm due to heating from the light emission part. Avoid long periods of contact between the distal end of the aScope 5 Cysto HD and the mucosal membrane, as sustained contact with the mucosal membrane may cause mucosal injury.
6. aScope 5 Cysto HD camera images must not be used as an independent diagnostic of any pathology. Doing so may result in incorrect or missing diagnosis. Physicians must interpret and substantiate any findings by other means and in the light of the patient's clinical characteristics.
7. Advance, withdraw or activate bending section of the aScope 5 Cysto HD with caution if an endoscopic instrument is protruding from the distal end of the working channel as this may cause mucosal injury.
8. Do not activate an energized endoscopic instrument (e.g. laser equipment, electrosurgical equipment) in the aScope 5 Cysto HD before the distal end of the instrument can be seen in the image on the displaying unit as this can lead to patient injury or damage aScope 5 Cysto HD.
9. Do not damage the insertion portion during use as it may leave parts of the product inside the patient or expose sharp surfaces that may cause damage to mucosa. Care should be taken to avoid damaging the insertion portion when using aScope 5 Cysto HD with endoscopic instruments.
10. Always watch the live image on the displaying unit when inserting or withdrawing the aScope 5 Cysto HD or operating the bending section. Looking at a recorded image may result in damage to mucosa or tissue.
11. Using electrosurgical equipment with aScope 5 Cysto HD may disturb image on displaying unit.
12. Do not use aScope 5 Cysto HD with laser equipment or electrosurgical equipment if flammable or explosive gases are present in the immediate area of aScope 5 Cysto HD as this can lead to patient injury, damage aScope 5 Cysto HD or disturb image on displaying unit.
13. Patient leakage currents may be additive and too high when using an energized endoscopic instrument in the aScope 5 Cysto HD. Only energized endoscopic instruments classified as "type CF" or "type BF" applied part shall be used with aScope 5 Cysto HD to minimize total patient leakage current. Failure to do so may cause major cardiac conductivity impairment and hemodynamic instability.
14. Irrigation by insufflation of air, inert gas prior to electrosurgery or using laser assist gas may cause a gas embolism leading to stroke or ischemia.
15. Do not use aScope 5 Cysto HD during defibrillation as this may result in electrical shock to the user.
16. When using compatible laser equipment, the user must be familiar with safety precautions, guidelines, and proper use of the laser equipment, including, but not limited to, proper eye and skin protection to avoid laser injuries.

## CAUTIONS

1. Have a suitable backup system readily available in case a malfunction should occur.
2. US federal law restricts this device for sale only by, or on the order of, a physician.
3. The aScope 5 Cysto HD may be inserted through a compatible access sheath. Advance gently to the desired treatment area and do not bend the controllable portion inside the access sheath.

4. Always progress with caution when inserting and advancing endoscopic instruments in a bent endoscope as this can cause damage to the working channel and introduce patient injury such as mucosal abrasion, bleeding and/or perforation.

## 1.7. Undesirable side effects

Undesirable side effects in relation to flexible cystoscopy (not exhaustive):

Intra-procedural pain or discomfort, hematuria, abdominal pain, dysuria – pain and discomfort on voiding, increased voiding frequency, urethral narrowing (strictures) due to scar tissue formation and urinary tract infections (UTI).

Potential adverse events in relation to antegrade flexible nephroscopy (not exhaustive):

mucosal abrasion, bleeding, perforation, sepsis, hematuria and urinary tract infections (UTI).

## 1.8. General notes

If, during the use of this device or as a result of its use, a serious incident occurs, please report it to the manufacturer and to your national authority.


## 2. Device description

aScope 5 Cysto HD is to be connected to an Ambu displaying unit. For information about Ambu displaying units, please refer to the Ambu displaying units' Instructions For Use.

### 2.1. Device parts

**Ambu® aScope™ 5 Cysto HD –  
Single-use device**

**Part numbers**

602001000 aScope 5 Cysto HD – Reverse Deflection  
603001000 aScope 5 Cysto HD – Standard Deflection

aScope 5 Cysto HD variants (#602001000 and #603001000) are not available in all countries. Please contact your local sales office.

### 2.2. Product compatibility

aScope 5 Cysto HD has been designed to be used in conjunction with:

#### Displaying units

- Ambu® aView™ 2 Advance
- Ambu® aBox™ 2

#### Endoscopic accessories and instruments

- Irrigation set (line and sterile water or saline bag) with Luer connection.
- Syringe and other Luer connecting accessories.
- Endoscopic instruments labelled for use in a working channel size of (ID) 2.0 mm/6.0 Fr or less\*.
- Access Sheath with an inner diameter of no less than 20 Fr.
- Holmium YAG laser (2.1 microns wavelength).
- Thulium fiber laser (1.92 – 1.96 microns wavelength).
- High frequency electrosurgical equipment fulfilling EN 60601-2-2. To keep high frequency leakage currents within allowed limits, the maximum sinus peak voltage level of the electrosurgical unit must not exceed 2.2 kVp.

\* There is no guarantee that instruments selected solely using this minimum working channel size will be compatible in combination. The compatibility of selected instruments should be tested before the procedure.


#### Lubricants and solutions

- Isotonic saline solution
- Iodine based (30 g) contrast agent suitable for cystoscopy and/or nephroscopy
- Water based soluble lubricants suitable for cystoscopy and/or nephroscopy

#### Other equipment

- X-ray

## 2.3. aScope 5 Cysto HD parts

No.	Part	Function
1	Handle	Suitable for left and right hand.
2	Control lever	Moves the distal end up or down in a single plane.
3	Working channel entry	Allows for instillation of fluids and insertion of endoscopic instruments.
	Working channel	Can be used for instillation of fluids and insertion of endoscopic instruments.
4	Insertion cord	Flexible insertion cord.
5	Bending section	Maneuverable part.
6	Distal end	Contains the camera, light source (two LEDs), as well as the working channel exit.
4-5-6	Insertion portion	The combination of the insertion cord, bending section, and distal end.
7	Connector on the aScope 5 Cysto HD cable	Connects to the grey socket on Ambu displaying units.
8	aScope 5 Cysto HD cable	Transmits the image signal to Ambu displaying units.
9	Protection pipe	Protects the insertion cord during transport and storage. Remove before use.

## 3. Use of aScope 5 Cysto HD

Numbers in gray circles below refer to illustrations on page 2.

### 3.1. Preparation and inspection of aScope 5 Cysto HD

#### Visual Inspection of the aScope 5 Cysto HD 1

1. Check that the pouch seal is intact before opening. Discard the aScope 5 Cysto HD if the pouch seal has been damaged, or the expiry date has been exceeded 1.
2. Make sure to remove the protection pipe from the insertion cord.
3. Check that there are no impurities on or damage to the aScope 5 Cysto HD such as rough surfaces, sharp edges or protrusions which may harm the patient 2.
4. Check the deflection of the controllable portion, by moving the control lever on the handle with your thumb to bend the bending section as much as possible. Then slide the control lever slowly to its neutral position. Confirm that the bending section functions smoothly and correctly and returns to a neutral position 5.

**Discard the aScope 5 Cysto HD if it fails any of the check points mentioned above. Refer to the Ambu displaying units' Instructions For Use for preparing and turning the Ambu displaying units on 3.**

### Inspection of the image

1. Turn the displaying unit on. Connect the aScope 5 Cysto HD to the Ambu displaying unit by plugging the connector on the aScope 5 Cysto HD cable marked with a grey arrow into the corresponding grey female connector on the Ambu displaying unit. Carefully align the arrows on the connector on the aScope 5 Cysto HD cable with the port on the Ambu displaying unit to prevent damage to the connectors **3 4**.
2. Verify that a clear and correctly oriented live video image appears on the Ambu displaying unit by pointing the distal end of the aScope 5 Cysto HD towards an object, e.g. the palm of your hand.
3. Adjust the image preferences on the Ambu displaying unit if necessary (please refer to the Ambu displaying unit's Instructions For Use).
4. If the object cannot be seen clearly, wipe the distal end of the aScope 5 Cysto HD using a sterile cloth.

### Preparation of aScope 5 Cysto HD

1. Test fluid instillation by connecting an infusion set or syringe with sterile water or saline solution with Luer connection directly to the working channel entry or via a stopcock. Ensure that there are no leaks and that water is emitted from the distal end **6**.

**Ensure that the check points above are successful before moving on to the procedure.**

## 3.2. Operating the aScope 5 Cysto HD


In the event of a malfunction during the cystoscopic or nephroscopic procedure, stop the procedure immediately, put the distal end of the aScope 5 Cysto HD in its neutral and non-angled position and slowly withdraw the cysto-nephroscope. Do not activate the control lever while withdrawing the cysto-nephroscope from the patient.

### Holding the aScope 5 Cysto HD and manipulating the distal end


The handle of the aScope 5 Cysto HD can be held in either hand. The hand that is not holding the cysto-nephroscope can be used to advance the insertion cord into the patient's lower urinary tract or kidney.

Use the thumb to move the bending lever. The bending lever is used to flex and extend the distal end of the cysto-nephroscope in the vertical plane.

- Depending on the bending lever, the model is called standard (lever up = tip up) or reverse (lever up = tip down).

Standard Deflection  
Lever Up = Tip Up  
Lever Down = Tip Down

Reverse Deflection  
Lever Up = Tip Down  
Lever Down = Tip Up

- The insertion cord should be held as straight as possible at all times in order to secure an optimal distal end bending angle.

### Insertion of the aScope 5 Cysto HD

For retrograde access lubricate the insertion cord with a soluble lubricant suitable for cystoscopy before the aScope 5 Cysto HD is inserted into the urethra. For antegrade access the aScope 5 Cysto HD may be inserted through a compatible access sheath. Gently advance the insertion cord according to standard practice for the patient anatomy. If the camera image of the aScope 5 Cysto HD becomes unclear, the distal end can be cleaned by withdrawing the cysto-nephroscope and cleaning the distal end.

### Aspiration and instillation of fluids 6

Aspiration may be required during the procedure. Prepare a syringe for this. When required, attach the syringe to the aScope 5 Cysto HD and apply an aspiration force to achieve the desired effect. For larger quantities of fluid, disconnect the syringe from the cysto- nephroscope, empty the syringe, and then reattach it to aspirate the remaining fluids.

Fluids e.g. sterile water or saline solution can be instilled through the working channel entry at the bottom of the aScope 5 Cysto HD handle by connecting a syringe or infusion set with Luer Lock connection directly to the working channel entry or via a stopcock. If using a sterile water or saline bag, make sure to place it so that potential spillage will not affect other equipment.

### Insertion of endoscopic instruments 7

If required, endoscopic instruments can be used with the aScope 5 Cysto HD.

Always make sure to select the correct size endoscopic instrument for the aScope 5 Cysto HD (see section 2.2.). Inspect the endoscopic instrument before using it. If there is any irregularity in its operation or external appearance, replace it. It is recommended not to have a fully deflected bending section when inserting the endoscopic instrument. Insert the endoscopic instrument into the working channel entry and advance it carefully through the working channel until it can be seen on the live image on the Ambu displaying unit. Do not activate endoscopic instruments inside the working channel. The distal end of the instrument should be seen in the image during use.

It should be recognized that the use of lasers may interfere with the normal endoscopic image and this interference is not necessarily indicative of a malfunction of the endoscopic system. A variety of factors can affect the quality of the endoscopic image during use of lasers. Factors such as intensity, high power setting, close distance of the instrument probe to the cysto-nephroscope end and excessive stone treatment can each adversely influence image quality.

### Withdrawal of the aScope 5 Cysto HD 8

When withdrawing the aScope 5 Cysto HD, make sure that the control lever is in the neutral position. Slowly withdraw the cysto-nephroscope while watching the live image on the displaying unit.

## 3.3. After use

### Visual check 9

Check for missing parts, evidence of damage, cuts, holes, sagging, or other irregularities on the bending section, distal end, or insertion cord of the aScope 5 Cysto HD. If required, take corrective action to determine whether any parts are missing and locate the missing part(s).

If corrective action is needed, follow the local hospital procedures. The elements of the insertion cord are visible in x-ray (radio opaque).

### Final steps

1. Disconnect the aScope 5 Cysto HD from the Ambu displaying unit.
2. Dispose of the aScope 5 Cysto HD, which is a single-use device 10. The aScope 5 Cysto HD is considered contaminated after use and must be disposed of in accordance with local guidelines for collection of infected medical devices with electronic components. The product design and materials used are not designed for reuse and cannot withstand the reprocessing procedures used for reprocessing endoscopes without risking degrading and contamination. Reprocessing/reuse could damage the scope and cause injury to the patient.

## 4. Technical product specifications

### 4.1. aScope 5 Cysto HD specifications

Insertion portion	aScope 5 Cysto HD	Optical system	aScope 5 Cysto HD
Bending angle <sup>1</sup>	210° ± 15° ↑ min. 120° ↓	Direction of view	0°

Insertion portion	aScope 5 Cysto HD	Optical system	aScope 5 Cysto HD
Insertion cord diameter	16.2 Fr $\pm$ 0.3 Fr/ 5.4 mm $\pm$ 0.1 mm (0.21" $\pm$ 0.004")	Field of view	120° $\pm$ 18°
Distal end diameter	Max 16.8 Fr Max 5.6 mm (Max 0.22")	Depth of field	3 – 100 mm
Maximum diameter of insertion portion	Max 18 Fr/ 6.0 mm (0.24")	Illumination method	LED
Working length	388 mm $\pm$ 12/-8 mm (15.3" $\pm$ 0.5"/-0.3")		
Working channel		Sterilisation	
Minimum working channel width <sup>2</sup>	Min. 6.6 Fr/ 2.2 mm (0.086")	Method of sterilisation	ETO
Storage and transportation condition		Operating environment	
Transportation temperature	-10 – 55 °C (14 – 131 °F)	Temperature	10 – 40 °C (50 – 104 °F)
Storage temperature <sup>3</sup>	10 – 25 °C (50 – 77 °F)	Relative humidity	30 – 85 %
Transportation relative humidity	10 – 95 %		
Storage relative humidity	10 – 85 %	Atmospheric pressure	80 – 106 kPa
Atmospheric pressure	50 – 106 kPa (7 – 15 psi)	Altitude	$\leq$ 2000 m

1. Please be aware that the bending angle can be affected if the insertion cord is not kept straight; or endoscopic instruments have been inserted.
2. There is no guarantee that accessories selected solely using this minimum working channel width will be compatible in combination.
3. Storage under higher temperatures may impact shelf life.




## 5. Troubleshooting

If problems occur with the system or product, please use this troubleshooting guide to identify the cause and correct the error.

Problem	Possible cause	Recommended action
No live image on the Ambu displaying unit, but the user interface is present on the Ambu displaying unit or the image shown is frozen.	The aScope 5 Cysto HD is not connected to the Ambu displaying unit.	Connect the aScope 5 Cysto HD to the grey port on the Ambu displaying unit.
	The Ambu displaying unit and the aScope 5 Cysto HD are not communicating correctly.	Restart the Ambu displaying unit (please refer to the Ambu displaying unit's Instructions For Use).
	The aScope 5 Cysto HD is damaged.	Replace the aScope 5 Cysto HD with a new one.
	A recorded image is shown.	Return to the live image (please refer to the Ambu displaying unit's Instructions For Use).

Problem	Possible cause	Recommended action
Low picture quality.	Unwanted fluids etc. on the distal end.	Gently rub the distal tip against the mucosa. If the lens cannot be cleaned this way remove the endoscope and wipe the lens with sterile gauze.
Absent or reduced flow of fluid, e.g. sterile water or saline solution, or difficulty inserting the endoscopic instrument through the working channel.	The working channel is blocked.	Clean the working channel using a cleaning brush or flush the working channel with sterile water or saline using a syringe. This should not be done with the aScope 5 Cysto HD inside the patient to avoid potential flush debris or other foreign matter into the patient.
	The bending section is not in neutral position.	Move the bending section into neutral position.

## 6. Symbol description

Symbols for the aScope 5 Cysto HD devices	Description	Symbols for the aScope 5 Cysto HD devices	Description
	Working length of the aScope 5 Cysto HD insertion cord		Single sterile barrier system
 Max OD	Maximum insertion portion width (maximum outer diameter)		Global Trade Item Number
 Min ID	Minimum instrument channel width (Minimum inner diameter)		Country of manufacture. Date of manufacture
	Field of view		UL mark on electronic products (UL Recognized Component Mark for Canada and the United States)
	Direct current input		UK Responsible Person
	Electrical safety type BF applied part		UK Conformity Assessed. Indicates that the product is in compliance with UK legislation for medical devices and has passed conformity assessment by UK Approved Body
	Medical device		
	Importer		

A full list of symbol explanations can also be found on [ambu.com/symbol-explanation](http://ambu.com/symbol-explanation).

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Ambu and other trademarks are trademarks of Ambu A/S
**Ambu A/S**
**Ambu Ltd**