


## CARDIOLOGY ACCESSORIES

### Features and Benefits

- Fixes electrode leadwires on patients during Holter and Stress test
- Reduces movement artefacts on the ECG
- Minimizes the risk of electrodes falling off
- MR safe
- Adhesive size: Ø 55 mm
- PVC free
- Latex free
- Biocompatible



### Skin Fix

Ref: 2703/50

Units/pouch	Units/inner box	Units/outer box
50	500	3000

### Features and Benefits

- Allows skin preparation before electrode placement
- Gently abrasive to reduce skin impedance
- Adhesive on one side, to use with one finger
- Small size fits finger tip
- Size: 19 x 28 mm
- PVC free



### Skin Prep

Ref: 2121M

Units/pouch	Units/inner box	Units/outer box
500	1500	15000

# CARDIOLOGY ACCESSORIES

## Features and Benefits

- Connects 3 mm pin cables to tab electrodes
- Ensures secure connection
- Easy to use
- Reusable
- PVC free



## Snap Clip

Ref: 0334M

Units/pouch	Units/inner box	Units/outer box
10	N/A	1000

## Features and Benefits

- Connects 4 mm banana pin to snap or tab electrodes
- Ensures secure connection
- Allows quick and easy one-hand application
- Reusable
- PVC free



## Multi Clip

Ref: 2952/10

Units/pouch	Units/inner box	Units/outer box
10	N/A	N/A

# Ambu

## Ambu Limited

Incubator 2,  
Alconbury Weald Enterprise Campus,  
Cambridgeshire,  
PE28 4XA

T 01480 498 403  
ambu.co.uk

LBL-006912 - V01 - 2022/07 - Ambu A/S. Technical data may be modified without further notice. TCC 11330

**CE** US: Rx only