

DATASHEET

The Ambu King Mask is an anatomically shaped, disposable anaesthetic facemask that is designed to support a tight seal around the facial anatomy with minimal pressure to the face. The Ambu® King Mask has a thin cushion which improves the ability to create a tight seal on a variety of patients. Furthermore, by inflating or deflating the cuff using the check valve in order to suit the patient's anatomical contours, the seal can be adjusted.

The Ambu King Mask is listed as MR Safe and is therefore suitable for use in an MR environment.

The Ambu® King Mask is also available in a scented version with Fresh Scent, Bubble Gum, Strawberry and Cherry scents. This can have a calming effect during induction, especially for children¹.

Five Scents

- Unscented (Blue)
- Fresh Scent (Green)
- Bubble Gum (Purple)
- Strawberry (Pink)
- Cherry (Red)

Seven Sizes

- Size 1
- Size 2
- Size 3
- Size 4
- Size 5
- Size 6
- Size 7RD

Ambu® King Mask

Disposable anaesthesia face mask



KEY BENEFITS

Seal can be enhanced by deflating or inflating the cushion via the check valve

Cushion is thin and pliable, for an excellent face seal

Ribbed on the edge of the cushion, for improved grip

Anatomically shaped cuff for improved adjustment to the face with minimal pressure

Clear dome for easy observation of the patient's condition








MR Safe

Available with or without hook ring*






¹ Fukumoto M, Arima H, Ito S, Takeuchi N, Nakano H. Distorted perception of smell by volatile agents facilitated inhalational induction of anaesthesia. Paediatr Anaesth. 2005 Feb; 15 (2): 98-101.

* Not available for sizes 1 and 2.

ORDERING SPECIFICATIONS WITHOUT HOOK RING

Mask	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7RD
Packaging Color Code							
Original (unscented)	1012UB	1022UB	1032UB	1045UB	1055UB	1065UB	1075UB
Fresh Scent (mint)	2012SB	2022SB	2032SB	2045SB	2055SB	2065SB	2075SB
Bubble Gum		1122SB	1132SB	1145SB			
Strawberry		1222SB	1232SB	1245SB			
Cherry		1322SB	1332SB	1345SB			
Scented masks variety pack		1423SB					

WITH HOOK RING

Mask	Size 3	Size 4	Size 5	Size 6	Size 7RD
Packaging Color Code					
Original (unscented)	1032UH	1045UH	1055UH	1065UH	1075UH
Fresh Scent (mint)	2032SH	2045SH	2055SH	2065SH	2075SH
Bubble Gum	1132SH	1145SH			
Strawberry	1232SH	1245SH			
Cherry	1332SH	1345SH			
Scented masks variety pack	1433SH	1443SH			

Part #'s ending in 2 = 20 per case. Part #'s ending in 5 = 50 per case
 Part #'s ending in 3 = 30 per case (Scented mask variety pack 10 of each of Bubble Gum, Strawberry, & Cherry)

MATERIALS

Product	Material
Cushion	PVC Polyvinyl Chloride (Phthalate-free)
Crown	Polycarbonate
Hook ring	Polycarbonate
Valve Body	Polycarbonate
Valve Stem	Silicone Rubber
Valve Plug	PP Polypropylene

Packaging

Pouch	LDPE
Box	Cardboard

Environment

Phthalate-free	Yes
Not made with natural rubber latex	Yes

SPECIFICATIONS

Size	Weight
1	15 g
2	18.0 g
3	33.0 g
4	39.0 g
5	48.0 g
6	59.0 g
7RD	53.0 g

Technical Information

Connector size	15 mm OD connector for #1 and #2 according to EN ISO 5356-1. 22 mm ID connector for #3, #4, #5, #6 and #7RD according to EN ISO 5356-1.
Check Valve	Compatible with Luer slip syringe
Operation temperature limits	-20 °C to +50 °C (-4 °F to +122 °F) according to EN 1789.
Storage temperature limits	-40 °C to +70 °C (-40 °F to +158 °F) according to EN 1789 and EN ISO 10651-4.

Recommended long term storage in closed packaging at room temperature, away from sunlight.



Unscented



Fresh Scent



Strawberry



Cherry



Bubble Gum

Ambu

Ambu Limited

Incubator 2
 Alconbury Weald
 Enterprise Campus
 Cambridgeshire
 PE28 4XA
 T +44 1480 498 403
 ambu.co.uk

CE 2797