

SPECIFICATIONS – AMBU REDI-ACE

Adjustable sizes	Neckless, Short, Regular and Tall
Number of settings	8 adjustments
Dimensions	554 x 215 x 14 mm (21.8" x 8.5" x 0.6")
Weight	200 g (7 ounces)
Color	Blue/light grey
Translucency	CT and MRI compatible
Storage temperature	-30° to +70° C
Operation temperature	-10° to +40° C



Ambu Redi-ACE Adult

SPECIFICATIONS – AMBU REDI-ACE MINI

Adjustable sizes	Infant, Pediatric and Small Adult
Number of settings	6 adjustments
Dimensions	447 x 188 x 14 mm (17.6" x 7.4" x 0.6")
Weight	100 g (3.5 ounces)
Color	Yellow/light grey
Translucency	CT and MRI compatible
Storage temperature	-30° to +70° C
Operation temperature	-10° to +40° C



Ambu Redi-ACE Mini

ORDERING INFORMATION

	CAT. NO
Redi-ACE Extrication Collar, 30/case	472 001 000
Redi-ACE Mini Extrication Collar, 30/Case	472 002 000



Flat storage



Locking buttons



Cannula hooks

Main benefits of Ambu Redi-ACE:

- **Easy to use**
Ambu Redi-ACE has an automatic flip up chin and two locking buttons. It can be placed on the patient using only 3 steps; sizing – locking – placement.
- **2 collars for all patients**
With 8 sizes of Ambu Redi-ACE and 6 sizes of Ambu Redi-ACE Mini the majority of patients can be treated.
- **Practical characteristics**
Ambu Redi-ACE can be stored completely flat. Additionally the Ambu Redi-ACE has cannula hooks to secure the oxygen tube.

Ambu® Redi-ACE Adult & Redi-ACE Mini

Description

The Ambu Redi-ACE is a unique adjustable collar. Combined with the Redi-ACE Mini these two collars offer 14 sizes ensuring that patients can have the most effective and accurate fit possible. Made from radio translucent materials, the collar is packaged completely flat and has a unique automatic flip up chin. Ambu Redi-ACE is a disposable adjustable collar for immobilizing the cervical column. The Ambu Redi-ACE incorporates latches that allow the rescuer to size the collar as accurate as possible to the victim's neck size. Two safety lock buttons hold the collar in place when adjusted.

MATERIALS – AMBU REDI-ACE

White/light grey plastic	Polypropylene
Foam	IXPE

MATERIALS – AMBU REDI-ACE MINI

Yellow plastic	Polypropylene
Foam	IXPE

Ambu A/S
Baltorpbakken 13
DK-2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
www.ambu.com

496 8500 01 - V01 - 06/2011 - Ambu A/S - Technical data may be modified without further notice.

Cervical immobilization using Ambu Redi-ACE*

Sitting position



Manual neutral in-line stabilization of the patient's head and neck, provided by first provider.



Second provider measures with the finger the patient's neck between the lower jaw and the shoulder.



Second provider adjusts the collar before measured size.



Second provider locks the collar in the proper size.



Second provider applies the collar while the first provider maintain the neutral in-line position of head and neck stabilisation.



Second provider secures the collar using the Velcro straps. Manual in-line stabilisation of the head and neck is maintained until the patient is secured to an immobilization device.



Supine position



* An extrication collar is only part of total patient immobilization. Maintain neutral head and neck alignment during the entire procedure. Ambu does not attempt to define medical practice. For detailed instructions on patient immobilization please see international and national guidelines and recommendation.