

Ambu A/S
Baltorpbakken 13
DK-2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 53
www.ambu.com

AMBU BUSINESS AREAS

- RESPIRATORY CARE
- CARDIOLOGY
- NEUROLOGY
- TRAINING
- IMMOBILIZATION









Product Catalogue

Ambu[∦]

- RESPIRATORY CARE
- CARDIOLOGY
- NEUROLOGY
- TRAINING
- IMMOBILIZATION











Ambu develops, produces and markets diagnostic and life-supporting equipment and solutions to hospitals and rescue services. It has five business areas: Respiratory Care, Cardiology, Neurology, Training and Immobilization, in which the most important products are ventilation products for artificial respiration and single-use electrodes for ECG tests and neurophysiological mappings.

Approx. 98% of the Group's products is sold on export markets through Ambu's foreign sales companies and our distributors.

The total number of employees is 1040 of whom 380 work in Denmark and 660 in countries outside Denmark.

Index

IMMOBILIZATION

Ambu Perfit	2
Ambu Mini Perfit ACE	
Ambu Perfit ACE	
Ambu Head Wedge	6
Ambu Military Perfit ACE	7
Amhu Military Head Wedge	-





Ambu[®] Perfit[™] Extrication Collars



THE PERFECT FIT FOR NEUTRAL ALIGNMENT

Ambu offers two variations of collar. The Ambu Perfit and the Ambu Perfit ACE (Adjustable Collar for Extrication). Both types of collar are designed to assist with the maintenance of neutral alignment, prevention of lateral sway and anterior-posterior flexion and extension of the cervical spine during transport and routine patient care or movement.

ALL AMBU EXTRICATION COLLARS HAVE THE FOLLOWING FEATURES

- Ambu collars are fully disposable, you no longer have to worry about collars being returned to the source of origin and cross contamination issues are avoided
- Ambu collars come with a pre-shaped flip-chin piece to ensure they can be stored completely flat. Furthermore the flip-chin piece facilitates intubation and removal of vomit from a patient's mouth

- Ambu collars come with a comprehensive sizing system for easy and accurate sizing
- One-piece extrication collars simplify application and at the same time there are no extra pieces that can get lost
- Nasal cannula holders are located on the collar body to secure the oxygen tubes (except for infant)
- All collars are radiotranslucent, CAT & MRI compatible and come with a ventilated posterior shell for fluid drainage and a large hole for easy access. All of these features allow for patient examinations to take place safely without removing the collar
- Safety buttons allow for maximum reliability



Ambu® Perfit

AMBU PERFIT EXTRICATION COLLAR IS A ONE-SIZE COLLAR,

which comes in four adult sizes, a pediatric size and an infant size.



TECHNICAL DATA

Shell m	naterial: polyethylene	Colour: White/light grey
Sizes:	1 (infant), 2 (paediatric), 3 (neckless), 4 (short),	
	5 (regular) and 6 (tall)	



ORDERING INFORMATION	CAT. NO.
Perfit Infant: size 1	264 201 000
Perfit Paediatric: size 2	264 202 000
Perfit Neckless : size 3	264 203 000
Perfit Short: size 4	264 204 000
Perfit Regular: size 5	264 205 000
Perfit Tall: size 6	264 206 000
Carry bag for collars, fits 6 pieces	264 008 000
Sample pack, one of each Perfit size	264 020 000

Ambu® Mini Perfit ACE™ Ambu® Perfit ACETM



28 SIZES IN 2 COLLARS

The Ambu Perfit ACE is a revolution in disposable adjustable collars. Combined with the Mini Perfit ACE, these two collars offer an incredible 28 sizes between them, ensuring that patients can now have the most effective and accurate fit ever.

The Ambu Perfit ACE is the only adjustable adult one-piece rigid cervical spine immobilisation

device with 16 precise settings for a customised fit. It adjusts within the 4 standard adult settings ranging from Neckless to Tall.

The Mini Perfit ACE is the only small adult/child adjustable one-piece rigid cervical spine immobilisation device with 12 precise settings for a customised fit. It adjusts from infant up to small adult neck sizes.



TECHNICAL DATA

Ambu Mini Perfit ACE		
Shell material: polyethylene	Number of adjustments: 12	
Adjustable sizes: infant, paediatric and small adult	Colour: Yellow/light grey	
Ambu Perfit ACE		
Shell material: polyethylene	Number of adjustments: 16	
Adjustable sizes: neckless, short, regular and tall	Colour: White/light grey	



ORDERING INFORMATION	CAT. NO.
Perfit ACE (1 pcs)	002 810 001
Mini Perfit ACE (1 pcs)	000 281 106

Ambu[®] **Head Wedge To be attached to the back board**







See interactive flash on our home page www.ambu.com

A REVOLUTION IN DESIGN

The unique Ambu Head Wedge is designed to immobilise a patient's head instantly and effectively when placed on a back board. It is applied easily and quickly in just four simple steps.

With its innovative construction the Head Wedge enables a perfect fit to the patient in one go, saving time at the scene when seconds count.

Made of plastic, it is superior to any cardboard product; there is no need to worry about bad weather because the Head Wedge is even efficient on wet back boards.

FEATURES

- Fully disposable
- Manufactured out of recyclable plastic, making it both waterproof and very durable
- It is totally radiolucent, CT and MRI compatible, ensuring that the patient is effectively immobilised even during examinations
- The Head Wedge stores completely flat saving valuable storage space, it can even be stored attached to the backboard
- Fully adjustable to fit all sizes from children to adults
- Fits all known backboards
- Strong, self-adhesive head strap with medical grade glue ensures a secure hold without any discomfort to the patient

ORDERING INFORMATION	CAT. NO.
Ambu Head Wedge (with two straps)	000 264 036
Ambu Head Wedge (with one strap)	000 264 034
Ambu Head Wedge (no strap)	000 264 035

Ambu[®] **Disposable Military Perfit ACETM**

A REVOLUTION IN DISPOSABLE ADJUSTABLE COLLARS

The Military Perfit ACE is a fully adjustable one piece cervical spine immobilisation device designed to give

a customised fit. It is identical to the standard Perfit ACE , but the Military perfit ACE has been designed in camoflauge colours to suit military specifications.



ORDERING INFORMATION

CAT. NO.

Ambu Military Head Wedge

000 281 000 MIL

Ambu® Military Head Wedge

INSTANT AND EFFECTIVE IMMOBILISATION

The fully adjustable, disposable Ambu Military Head Wedge immobilises a patients head instantly and

effectively in four simple steps. It is similar in design to the standard Ambu Head Wedge but it has been manufactured in camoflauge colours as per military specifications.



ORDERING INFORMATION

CAT. NO.

Ambu Military Head Wedge

000 264 064