

# SINGLE-USE EFFICIENCY, NEXT-LEVEL IMAGING

Ambu® single-use cystoscopy solutions



**Ambu**



**What if you could combine the efficiency of your Ambu® aScope™ 4 Cysto solution with a superior HD image?**

**Now you can with the new aScope™ 5 Cysto HD.**

**Introducing Ambu® aScope™ 5 Cysto HD**

aScope 5 Cysto HD gives you the efficiency you're familiar with from your aScope 4 Cysto endoscope, but with a superior HD image.

Just as with aScope 4 Cysto, it is single-use, so there is no deterioration affecting image quality, and you don't have to worry about damaging your aScope 5 Cysto HD endoscope.

**High image quality using HD technology**

The combination of a high-resolution camera chip in aScope 5 Cysto HD and our most innovative full-HD display unit, Ambu® aView™ 2 Advance, delivers a superior image.

With aScope 5 Cysto HD and the fully compatible aView 2 Advance, you can count on clear sharp images that make it easy to identify anatomical structures for a wide range of procedures in outpatient clinics, patient wards, ICUs and emergency rooms. For operating room procedures, aScope 5 Cysto HD is fully compatible with Ambu® aBox™ 2 as well.

Use the table on the next page as a guide to how aScope 5 Cysto HD can best support your current requirements.

**WANT TO LEARN MORE AND TRY  
Ambu® aScope™ 5 Cysto HD?**

Contact us with for further details or to arrange a trial:  
[ambu@ambu.com](mailto:ambu@ambu.com)

Or use the guide on the next page to see how  
aScope 5 Cysto HD matches your clinical needs

# aScope 5 Cysto HD

# HOW THE SOLUTION MATCHES YOUR CLINICAL NEEDS

aScope 4 Cysto	WHEN YOUR WORKPLACE NEEDS	aScope 5 Cysto HD
<p>✓</p> <p>✗</p>	<p><b>IDEAL APPLICATION</b></p> <p>A cost-efficient solution for in patient consults, the emergency room and the outpatient clinic</p> <p>Superior image quality for the operating theatre, including Advanced Red Contrast (ARC) enhancement algorithm</p>	<p>✓</p> <p>✓</p>
<p>✓</p> <p>✗</p> <p>✗</p>	<p><b>PROCEDURES</b></p> <p>Therapeutic cystoscopies (e.g. ureteral stent removal, bladder stone treatment and biopsies)</p> <p>Diagnostic cystoscopies (e.g. follow-up cystoscopies)</p> <p>Complex cystoscopies (e.g. treatment with HF or laser tools)</p>	<p>✗</p> <p>✓</p> <p>✓</p>
<p>✓</p> <p>✗</p>	<p><b>DISPLAY UNIT COMPATIBILITY</b></p> <p>Ambu® aView™ 2 Advance</p> <p>Ambu® aBox™ 2</p>	<p>✓</p> <p>✓</p>
<p>✓</p> <p>✗</p>	<p><b>DISPLAY UNIT PLACEMENT</b></p> <p>Set it up on a flat surface using the built-in stand or mounted on an IV pole</p> <p>Fits neatly in the tower in your operating theatre</p>	<p>✓</p> <p>✓</p>
<p>✓</p> <p>✗</p>	<p><b>DISPLAY UNIT ADJUSTMENT</b></p> <p>Rotate 180-degrees to connect in the best position for the given clinical situation</p> <p>Multi-directional: Can be rotated, raised or lowered to meet your viewing needs during the procedure</p>	<p>✓</p> <p>✓</p>
<p>✓</p> <p>✓</p>	<p><b>DOCUMENTATION VIA DISPLAY UNIT</b></p> <p>Ambu monitor connects via DICOM and PACS to enable you to seamlessly transfer your recorded patient imaging data</p> <p>Connects directly to your EMR documentation system to enable access to imaging data from the hospital's electronic medical records</p>	<p>✓</p> <p>✓</p>

## **FORWARD-THINKING SOLUTIONS**

Ambu has built on more than a decade of single-use endoscopy know-how to place a reliable, sterile and cost-effective solution in your hands. We launched the world's first single-use flexible bronchoscope in 2009, and in 2020, over 1 million of our single-use endoscopes were used in more than 6,000 hospitals globally.

### **The largest single-use R&D team in the world**

Ambu has hundreds of engineers in five innovation centres around the world - three of which focus exclusively on endoscopy solutions. They collaborate with medical specialists and personnel at our advanced production facilities. Our modular design approach, infrastructure and vast experience enable us to develop high-quality single-use technology across clinical areas and to deliver it efficiently.

### **Fully focused on single-use endoscopy**

Twenty per cent of our single-use revenue is fed back into product innovation, and in the next three years we are poised to launch dozens of single-use endoscopes across a range of clinical areas, including solutions for urology. With one experienced, reliable partner, you can standardise your single-use endoscopy solution across multiple clinical speciality areas.

**For more information, please visit [ambu.co.uk/urology](https://www.ambu.co.uk/urology)**

Ambu Limited  
Incubator 2  
Alconbury Weald  
Enterprise Campus  
Cambridgeshire  
PE28 4XA  
T +44 1480 498 403  
[ambu.co.uk](https://www.ambu.co.uk)