



Ambu®

Airway

# Ambu

Ambu's company philosophy has always been the desire to make life easier for doctors and to find innovative products that can save lives. Innovation and entrepreneurship are key for Ambu. Always have been. Always will be.



## Ambu 75 years – past, present and future

Ideas may be conceived in one person's head, however, it takes a collaboration of minds to see them evolve into real products.

Since 1937, we have been leading the way in healthcare. We develop ideas whether they arrive through meticulous research or by a clever question asked by a customer. All that matters is that we transform ideas into efficient solutions that work to our customer's advantage and improve patient safety.

### We listen to your needs

Our desire to listen to the needs of healthcare professionals has been the driving force behind the development of the Ambu bag, the Blue Sensor® electrodes or the world's first single-use videoscope, the aScope™.

Ambu devices are used every day to rescue, revive, diagnose, monitor or train. We are proud of this. But our products do nothing on their own. They are only useful in the skilled hands of doctors, nurses and paramedics around the world.

### Next 75 years

Our customers tell us that they need reliable high-quality products ready for immediate use. They expect excellent service and on-time delivery. They need innovation to help improve patient outcomes and optimise efficiency. They want training and clinically meaningful evidence.

Above all, our customers want us to truly understand the challenges of their patient care and to collaborate on the development of future solutions.

We are determined to continue to develop best-in-class products and the services you expect for your work each and every day. We do it now and we will do so for many decades to come.

# Fields of Excellence

## Anaesthesia

With the innovative development of the first self-inflating resuscitator, the "Ambu bag" in 1956, Ambu forever changed the history of saving lives. Now more than half a century later we are still dedicated to the development of innovative products that help save lives and simplify workflow and the daily life of the anaesthetists to the benefit of the patients.

Within the field of anaesthesia Ambu offers a wide range of products from resuscitators, face masks and laryngeal masks to the highly innovative single use flexible intubation scope. Products that all have their own place in the difficult airway algorithm.

## Emergency Care

We consistently focus on developing new devices to improve patient care and make the work of our customers easier. In the field of emergency care we offer a broad assortment of products, from resuscitators, suction pumps and CPR barriers to our well-known training manikins and immobilisation devices.

## Patient Monitoring & Diagnostics

In 1971, a Danish engineer by the name of Jan Rasmussen recognized the need to obtain reliable traces from biosignals. Since then we have been leading the way in the design and development of high quality single-use electrodes for accurate diagnosis.

Ambu offers a complete range of electrodes for different applications to meet your needs.

## Ambu's commitment to quality

Ambu products are designed and manufactured to comply with all relevant internationally recognised standards. Our entire organisation is committed to providing products that match the high expectations of users and patients, fully meet the requirements defined by standards and guidelines.

The Ambu quality management system operates under ongoing consideration for maintaining ISO 13485 quality management system certification, and is subject to regular inspections by the regulatory authorities in compliance with the international standard for medical device manufacturers.

Our commitment to quality includes an ongoing effort to improve on feedback, and forward product requests from our customers.



# Laryngeal Masks

The Ambu Aura series of laryngeal masks is composed of:

## AuraGain

Ambu AuraGain is the only anatomically curved SGA with integrated gastric access and intubation capability, taking patient safety and airway management efficiency to a new level.

## Ambu AuraOnce

Available in eight convenient sizes, AuraOnce from Ambu is the ideal alternative to a traditional face mask for achieving and maintaining control of the airway during routine and emergency anaesthetic procedures in fasted patients. It is also an excellent option when unexpected difficulties arise in connection with airway management.



## Features

### Easy to install "Easy glide"

Studies have concluded that the Ambu Aura Once laryngeal masks are easier to insert than other laryngeal masks. The special curvature of the tube is a faithful copy of the natural human anatomy. The user is not required to place his fingers in the mouth of the patient.

### Perfect seal

Clinical studies have shown that the Ambu laryngeal masks have a higher oropharyngeal seal. Specific studies indicate that the pressure of the cuff is not affected by the position of the head of the patient. Good sealing is crucial for ventilation, so if the pressure cannot be maintained, the ventilation of the patient will fail.

## Ambu AuraFlex

Ambu AuraFlex is a disposable, flexible laryngeal mask which is specially designed for ENT, ophthalmic, dental and other head and neck surgeries.

## Ambu AuraStraight

Ambu AuraStraight provides users with an even greater choice and an alternative to the already popular AuraOnce curved laryngeal mask for use both in anaesthesia and emergency medicine. Ambu AuraStraight comes in a full range of 8 sizes from neonate to large adult.

### Minimal trauma

Ambu laryngeal masks are less traumatic with regard to complications during anaesthesia pre-peri and postoperative procedures. Clinical studies conclude that postoperative complaints are rare.

### Versatile

Met with a wide range of products on offer, there is an Ambu laryngeal mask for virtually every need in all categories. Ambu develops and produces in close cooperation with healthcare professionals to provide the highest performance with innovative features.

Visit [www.ambu.co.uk](http://www.ambu.co.uk) for clinical studies.

## Summary of various laryngeal masks



### AuraGain

- DEHP free
- 8 sizes (1-6) for patients <5 to > 100kg
- Broad clinical use. Gastric access channel for managing gastric content
- Sterile with 3 year shelf life
- MR safe



### AuraOnce

- DEHP free
- 8 sizes (1-6) for patients <5 to > 100kg
- Broad clinical use
- Sterile with 3 year shelf life
- MR safe



### AuraFlex

- DEHP free
- 6 sizes (2-6), patients 10-20kg to > 100kg
- Great for dental surgery, tonsillectomy, or other oral or facial procedures where the airway tube should be moved to gain access
- Sterile with 3 year shelf life
- Not MR safe due to wire steel to the tube to strengthen



### AuraStraight

- DEHP free
- 8 sizes (1-6), patients <5kg to > 100kg
- Great for patients with a long neck, receding lower, or in patients lateral position
- Sterile with 3 year shelf life
- MR safe

Ambu® AuraGain™



Ambu AuraGain is the only anatomically curved SGA with integrated gastric access and intubation capability, taking patient safety and airway management efficiency to a new level.

The AuraGain is Ambu's 2nd generation laryngeal mask, satisfying 3 fundamental airway management needs by integrating gastric access and intubation capability in an anatomically curved single-use device that facilitates rapid establishment of a safe airway.

Features and Benefits:

- Integrated gastric access channel for managing gastric content
- The original anatomical curve, ensuring rapid placement
- Intubating capability using standard ET-tubes
- Integrated bite absorption area prevents airway occlusion
- Pilot balloon identifies mask size and provides tactile indication of degree of inflation
- Navigation marks for guiding flexible scope
- MR safe
- Phthalate-free material. Does not contain natural rubber latex

AuraGain  
Product Specification

	Infant		Paediatric		Adult			
Product item no.	408100000	408150000	408200000	408250000	408300000	408400000	408500000	408600000
Mask size	#1	#1½	#2	#2½	#3	#4	#5	#6
Packaging colour code								
Patient weight	< 5 kg	5-10 kg	10-20 kg	20-30 kg	30-50 kg	50-70 kg	70-100 kg	>100 kg

Ambu® AuraOnce™



Available in eight convenient sizes, AuraOnce from Ambu is the ideal alternative to a traditional laryngeal masks for achieving and maintaining control of the airway during routine and emergency anaesthetic procedures in fasted patients. It is also an excellent option when unexpected difficulties arise in connection with airway management.

The Ambu AuraOnce features a special 70° curve that carefully replicates natural human anatomy. Moulded directly into the tube for additional safety, the unique

shape of the mask makes correct insertion fast and easy without putting extra stress on the upper jaw. Moreover, like all Ambu laryngeal masks, AuraOnce has a higher oropharyngeal seal pressure to ensure reliable ventilation.

AuraOnce has also proven to be particularly useful in establishing a clear airway during resuscitation in profoundly unconscious patients with absent glossopharyngeal and laryngeal reflexes who may need artificial ventilation.

Features and Benefits:

- Ergonomically shaped airway tube for firm and ergonomical grip during insertion
- Soft cuff ensures a good seal with low mucosal pressure
- Reinforced tip resists folding during insertion and plugs the upper esophageal sphincter
- Cuff and airway tube moulded as single unit for extra safety – no separation
- Pilot balloon identifies mask size and provides precise tactile indication of degree of inflation
- Convenient depth marks for monitoring correct position
- Packaged sterile and ready for use in colour-coded pouches
- MR safe
- Phthalate-free material. DEHP free. Does not contain natural rubber latex

AuraOnce  
Product Specification

	Infant		Paediatric		Adult			
Product item no.	321100000	321150000	321200000	321250000	321300000	321400000	321500000	321600000
Mask size	#1	#1½	#2	#2½	#3	#4	#5	#6
Packaging colour code								
Patient weight	< 5 kg	5-10 kg	10-20 kg	20-30 kg	30-50 kg	50-70 kg	70-100 kg	>100 kg



Ambu® AuraStraight™



The Ambu® AuraStraight™ provides users with an even greater choice and an alternative to the already popular AuraOnce curved laryngeal mask for use both in anaesthesia and emergency medicine. Ambu AuraStraight comes in a full range of 8 sizes from neonate to large adult.

Features and Benefits:

- Reinforced tip resists folding over during insertion and plugs the upper esophageal sphincter
- Cuff and airway tube moulded as single unit for extra safety - no separation
- The surface has an EasyGlide texture to ease insertion
- Extra soft cuff ensures the best possible seal with least possible mucosal pressure
- Less invasive than endotracheal tube
- Pilot balloon identifies mask size and provides precise tactile indication degree of inflation
- Convenient depth marks for monitoring correct position of the mask
- Packaged sterile and ready for use
- Colour coded pouch with directions for use
- Available in 8 sizes fitting a wide range of patients
- MR safe and compatible
- Phthalate-free material. DEHP free. Does not contain natural rubber latex

AuraStraight  
Product Specification

	Infant		Paediatric		Adult			
Product item no.	324100000	324150000	324200000	324250000	324300000	324400000	324500000	324600000
Mask size	#1	#1½	#2	#2½	#3	#4	#5	#6
Packaging colour code	Orange	Grey	Red	Blue	Green	Yellow	Dark Blue	Red
Patient weight	< 5 kg	5-10 kg	10-20 kg	20-30 kg	30-50 kg	50-70 kg	70-100 kg	>100 kg

Ambu® AuraFlex™



Ambu® AuraFlex™ is a disposable, flexible laryngeal mask which is specially designed for ENT, ophthalmic, dental and other head and neck surgeries.

The challenge of a flexible laryngeal mask is the combination of needs: the ease of insertion and security of the airway needed by the anaesthesiologist, the

flexibility and room to work needed by the surgeon and the need to protect the patient from residual proteins. Ambu® AuraFlex™ combines these needs in an innovative, single-use solution.

Features and Benefits:

- High flexibility of the reinforced airway tube enables positioning away from the surgical field without loss of seal
- Integrated pilot tube eases the surgeons work and reduces the risk of cutting the pilot tube during surgery
- Kink-free wire-reinforced airway tube eliminates the risk of airway tube occlusion
- Full range from size 2 through size 6 - widest range of possible uses
- Thin and soft feel of the cuff with EasyGlide surface for easy and smooth insertion and removal
- Reinforced tip eases insertion and prevents folds during insertion - easy, fast and accurate positioning
- Depth marks for visual control of correct positioning of the mask
- Packaged sterile and ready for use
- Colour-coded packaging with directions for use
- Phthalate-free material. DEHP free. Does not contain natural rubber latex

AuraFlex  
Product Specification

	Paediatric			Adult		
Product item no.	327200000	327250000	327300000	327400000	327500000	327600000
Mask size	#2	#2½	#3	#4	#5	#6
Packaging colour code	Red	Blue	Green	Yellow	Dark Blue	Red
Patient weight	10-20 kg	20-30 kg	30-50 kg	50-70 kg	70-100 kg	>100 kg

# Ambu® Cuff Pressure Gauge



The Ambu® Cuff Pressure Gauge is designed to inflate and monitor the pressure of laryngeal masks and intubation tubes with low-pressure cuffs.

The dial has two zones of measurement graded in cm H<sub>2</sub>O with colour coding for optimum visibility. The Luer connection allows the gauge to be connected to laryngeal masks or intubation tubes, either directly or via the connecting tube.

The Ambu Cuff Pressure Gauge has an ergonomic design and is easy to operate thanks to the integrated inflation bulb and release valve. It addresses the particular needs of anaesthesia and resuscitation services.

## Features and Benefits:

- Simple and robust design
- Compatible with all Ambu laryngeal masks
- H<sub>2</sub>O graded in cm
- Easy-to-use inflation bulb and release valve
- Large dial for optimal display
- Standard Luer connector
- Phthalate-free material. DEHP free. Does not contain natural rubber latex

### Cuff Pressure Gauge Product Specification

Product item no.	321900000	321000702
Type	Ambu Cuff Pressure Gauge for inflating and testing cuffs	Tube for Cuff Pressure Gauge, 1 m

# Face Masks

## Ambu Face Masks

Ambu has a long history in the production and selling of quality masks for use both within the hospital and pre-hospital. All Ambu face masks are compatible with the different models and sizes of our resuscitators.

There are several options for choosing the correct facial mask according to the needs of the user and the patient. We have a fully transparent, reusable and disposable version of the classic blue silicone masks. The masks are typically used to support ventilation during the preparation of anaesthesia of the patient for the operation.

The Ambu face masks meet all-important features that users need, providing an airtight seal with minimal applied pressure. The mask is transparent for visual inspection of bleeding, vomiting and spontaneous breathing. The design and the shape of the cap is assured to support user manipulation during long-term ventilation. All our masks conform to the ISO rules relation to the size of the connections, and are compatible with all models of ventilation balloons on the market.

They have always been a great addition to the already existing Ambu resuscitators. Quality, effectiveness and comfort is recognised by all users. Ambu has expanded with the acquisition of King Systems.

Ambu King Masks are designed for comfort and safety, and are "The leading mask design". The edge is extremely flexible and feels sticky for excellent seal on the patient. These masks are superior in all shapes and sizes.





Summary of face masks



Ambu® UltraSeal disposable face mask

- DEHP free
- 7 sizes
- Specifically designed for use within anaesthetic departments
- Disposable
- With or without control valve
- 3 year shelf life



Ambu® King Mask disposable face mask

- Market Leading Design
- Cushion is pliable and tacky to the touch, for an excellent face seal
- Ribbed on top of the cushion, for a no-slip grip
- Anatomically shaped cuff enabling a tight seal with minimum applied pressure
- Flexible dome
- Clear dome for easy observation of the patient's condition
- Latex free
- MR Safe
- Phthalate-free material
- Color-coded hook rings for easy size identification\*



Ambu® Disposable face mask

- DEHP free
- 6 sizes
- These disposable masks are designed for use with both manual and automatic resuscitators
- Disposable
- With or without control valve
- 3 year shelf life



Ambu® Disposable open cuff face mask

- Transparent for visual check of bleeding, vomitus, and spontaneous breathing
- Available in 7 different colorcoded sizes for easy size identification
- Soft cuffs and shaped for a comfortable fit for both adults and children
- The masks are phthalate, pvc free, and are not made with natural rubber latex

\* Not available for sizes 1 and 2.

Ambu® UltraSeal face masks



The Ambu UltraSeal face mask is an anatomically correct disposable face mask designed specifically for anaesthetic departments.








It is a single use disposable face mask that comes both with and without a check valve. Sizes 3, 4, 5 and 6 come with a hook ring for attachment of the mask with a strap. The current range consists of 4 paediatric sizes, and 3 adult sizes up to size 6 Large Adult.








Features and Benefits:

- Transparent and crystal clear dome secures easy observation of the patient's condition during use
- Very soft anatomically shaped cuff that enables a tight seal with a minimum of pressure applied.
- Colour code on the package for quick and easy size identification
- All sizes are individually delivered in a clean medical box of 20 pieces total , or as single units.
- Hook ring on all versions apart of the neonate and infant dimensions. The hook ring can simply be removed if not needed
- Phthalate-free material. DEHP free. Does not contain natural rubber latex

Ambu® UltraSeal face masks

Product Specification

Ambu® UltraSeal face masks		Paediatric				Child/ Small Adult	Medium Adult	Large Adult
Product item no.	305102000	305104000	305106000	305108000	305010000	305012000	305014000	
Mask size	#0	#1	#2	#3	#4	#5	#6	
Packaging colour code								

Ambu® UltraSeal face masks with checkvalve		Paediatric				Child/ Small Adult	Medium Adult	Large Adult
Product item no.	305101000	305103000	305105000	305107000	305011000	305013000	305015000	
Mask size	#0	#1	#2	#3	#4	#5	#6	
Packaging colour code								

Ambu® King Mask disposable face masks



Because of the innovative quality design of the Ambu® King Mask, it is the market leading anaesthesia mask design available today. The Ambu® King Mask has a thin and tacky cushion which improves the ability to create a ventilation seal on a variety of patients. The design with rounded edges secures a smooth grip and user comfort. The Ambu® King Mask is made of phthalate-free PVC

material for patient and user safety. Furthermore, the Ambu® King Mask can be used in a MR environment, as it is MR Safe.

The Ambu® King Mask is also available in a scented version with Fresh Scent, Bubble Gum, Strawberry and Cherry scents. The Sweet Dreams™ masks release soothing aroma from natural oils in mask cushion.

Features and Benefits:

- Cushion is extremely pliable and tacky to the touch, for an excellent face seal
- Anatomically shaped cuff enabling a tight seal with minimum applied pressure
- Clear dome for easy observation of the patient's condition
- Ribbed on top of the cushion, for a no-slip grip
- Flexible, self-expanding dome
- Color-coded hook rings for easy size identification\*

\*Not available for sizes 1 and 2

Ambu® King Mask disposable face masks

Product Specification

Ambu® King Mask disposable face masks	Neonate	Infant	Toddler	Child/ Small Adult	Medium Adult	Large Adult	XL Adult
Mask size	#1	#2	#3	#4	#5	#6	#7RD
Packaging colour code	■	■	■	■	■	■	■
Original (unscented)	1012	1022	1032	1045	1055	1065	1075
Fresh Scent (mint)	2012	2022	2032	2045	2055	2065	2075
Bubble Gum		K1122	1132	1145			
Strawberry		1222	1232	1245			
Cherry		1322	1332	1345			
Scented masks variety pack		1423	1433	1443			

Part #'s ending in 2 = 20 per case • Part #'s ending in 5 = 50 per case  
Part #'s ending in 3 = 30 per case (Scented mask variety pack 10 of each of Bubble Gum, Strawberry, & Cherry)

Ambu® Disposable face masks



The Ambu disposable face mask is an anatomically shaped mask designed for use with manual and automatic fans and ventilators.

The mask comes in an assortment of 6 different sizes, and comes with or without a control valve, providing a good seal and minimum pressure on the patient.

Features and Benefits:

- Thumb Support for user convenience
- Soft, moulded, air filled edge to tightly fit on the face of the patient
- Crystal clear dome for easy observation of the condition of the patient
- For use with a single patient. Time-consuming cleaning, sterilisation and disinfection procedures are not necessary. The potential risk of cross infection is excluded.
- Hook ring on all versions apart from the neonate and infant dimensions. The hook ring can be removed if not needed
- Phthalate-free material. DEHP free. Does not contain natural rubber latex

Ambu® Disposable face masks

Product Specification

Ambu® Disposable face masks	Neonate	Infant	Toddler	Child/ Small Adult	Medium Adult	Large Adult
Product item no.	000252051	000252052	000252053	000252054	000252055	000252056
Ambu® Disposable face masks with checkvalve	Neonate	Infant	Toddler	Child/ Small Adult	Medium Adult	Large Adult
Product item no.	000252951	000252952	000252953	000252954	000252955	000252956



# Ambu® Disposable open cuff face masks



Ambu® Disposable Open Cuff Face Masks are designed for use with manual and automatic resuscitators and ventilators. Transparent domes for visual check of bleeding, vomitus, and spontaneous breathing.

The masks are made with soft cuffs and shaped for a comfortable fit for both adults and children. The masks are phthalate, pvc free, and are not made with natural rubber latex.

## Features and Benefits:

- Transparent for visual check of bleeding, vomitus, and spontaneous breathing
- Available in 7 different colorcoded sizes for easy size identification
- Soft cuffs and shaped for a comfortable fit for both adults and children
- The masks are phthalate, pvc free, and are not made with natural rubber latex

## Ambu® Disposable open cuff face masks

### Product Specification

Ambu® Disposable open cuff face masks	Neonate	Infant	Paediatric	Child/ Small Adult	Medium Adult	Large Adult	Extra Large Adult
Product item no.	306100000	306200000	306200000	306400000	306500000	306600000	306700000
Mask size	#1	#2	#3	#4	#5	#6	#7
Packaging colour code	■	■	□	■	■	■	■

# Additional information:

The Ambu and King Systems face masks are 100% latex and DEHP free.  
Most of our masks are available with and without a check valve.



# Resuscitators

With the innovative development of the new self-inflating "Ambu bag" in 1956, Ambu has forever changed the history of aid. Today, after more than half a century, Ambu is still involved in the development of innovative products that help save lives and simplify working conditions for anaesthetists and the benefit of patients.

In the field of ventilation Ambu offers a wide range of products, from the ventilator, face masks and laryngeal masks to flexible disposable videoscopes and video laryngoscopes. These products all have their place in the difficult airway algorithm.

## Summary of Ambu's disposable resuscitator



### SPUR II disposable resuscitator

- Available for adults, paediatrics and Newborn
- Material SEBS
- Disposable
- Not jacketed



## Ambu® SPUR II disposable resuscitator



The invention of the "Ambu bag" resuscitator in 1956 started Ambu on its continuous journey developing products that others may not consider possible.

Today the term "Ambu bag" is used all over the world as a generic term for a resuscitator, however only Ambu manufacture the real "Ambu bag"

The Ambu SPUR II (Single Patient Use Resuscitator) is designed for manual ventilation of neonates through to adults and now with M-Port for direct capnography recording.

### Features and Benefits:

- DIRECT etCO<sub>2</sub> Capability
- Integrated handle for reliable functionality
- Features our SafeGrip™ surface to ensure a secure grip onto the bag
- Unique shape is small and easy to hold. The ergonomic, lightweight design of Ambu SPUR II is made for optimal user handling and support even during extended ventilations
- Phthalate-free material. DEHP free. Does not contain natural rubber latex

### Ambu® SPUR II disposable resuscitator

#### Product Specification

	Adult	Paediatric	Infant
Stoke volume	one hand 800 ml/ two hands 1100 ml	450 ml	150 ml
Resuscitator volume	1475 m	635 ml	220 ml
Dimensions (length x diameter)	295 x 127 mm	234 x 99 mm	168 x 71 mm
Pressure-limiting valve	Pressure limiting valve*: 4.0 kPa (40 cm H <sub>2</sub> O)	4.0 kPa (40 cm H <sub>2</sub> O)	4.0 kPa (40 cm H <sub>2</sub> O)
MediPort	Standard Luer LS 6		

\*Adult version also available without pressure limiting valve





Adult	Paediatric	Infant
Item no: 325 002 000 Ambu SPUR II Adult, with patient valve, O <sub>2</sub> bag reservoir and disposable face mask size 5 (12 pcs. per box)	Item no: 330 002 000 Ambu SPUR II Paediatric, with pressure limiting valve, patient valve, O <sub>2</sub> bag reservoir and disposable face mask size 0 (12 pcs. per box)	Item no: 335 002 000 RH Ambu SPUR II Infant with pressure limiting valve, patient valve, O <sub>2</sub> tube reservoir and disposable face mask size 0 (12 pcs. per box) <b>with expiry date</b>
Item no: 325 014 000 Ambu SPUR II Adult, with pressure limiting valve, patient valve, O <sub>2</sub> bag reservoir and disposable face mask 1 x adult size and 1 x toddler size (12 pcs. per box)	Item no: 330 006 002 Ambu SPUR II Paediatric, with pressure limiting valve, O <sub>2</sub> bag reservoir and disposable face mask size 1 and size 2 (6 pcs. per box) <b>with expiry date</b>	Item no: 335 003 000 Ambu SPUR II Infant with pressure limiting valve, patient valve, O <sub>2</sub> tube reservoir and disposable face mask size 1 (12 pcs. per box)
Item no: 325 002 002 Ambu SPUR II Adult, with patient valve, O <sub>2</sub> bag reservoir and disposable face mask size 5 (12 pcs. per box) <b>with expiry date</b>	Item no: 330 006 000 Ambu SPUR II Paediatric, with pressure limiting valve, O <sub>2</sub> bag reservoir and disposable face mask size 1 and size 2 (6 pcs. per box, expanded)	Item no: 335 102 000 Ambu SPUR II Infant with pressure limiting valve, patient valve, O <sub>2</sub> bag reservoir and disposable face mask size 0 (12 pcs. per box)
Item no: 325 003 002 Ambu SPUR II Adult, with patient valve, O <sub>2</sub> bag reservoir and disposable face mask size 5 (6 pcs. per box) <b>with expiry date</b>		Item no: 335 002 000 Ambu SPUR II Infant with pressure limiting valve, patient valve, O <sub>2</sub> tube reservoir and disposable face mask size 0 (12 pcs. per box)
Item no: 325 012 002 Ambu SPUR II Adult, with pressure limiting valve, patient valve, O <sub>2</sub> bag reservoir and disposable face mask size 5 (12 pcs. per box) <b>with expiry date</b>		

# Ambu® PEEP Valves - Disposable

Ambu® PEEP Valves are designed for use with manual resuscitators or ventilators, where specified by the manufacturer.



## Features and Benefits:

- MR conditional, up to 3 Tesla (only disposable PEEP valve)
- Unique resistance
- Direct connection without adapter

## Ambu® PEEP Valves - Disposable

### Product Specification

Product item no.	199 103 001
Type	Ambu Disposable PEEP Valve 20, with adapter for respirators/CPAP system 22/19 mm

Can be purchased singly or in a box quantity of 10.

# Ambu® Disposable Pressure Manometer

The Ambu Disposable Pressure Manometer for single patient use is intended to be used for monitoring the patient's airway pressure.



## Features and Benefits:

- The slim design makes it easy to monitor pressure with minimum obstruction of the view of the patient's face
- As a single patient use product, the Ambu Disposable Pressure Manometer minimises the risk of cross infections
- Phthalate-free material. DEHP free. Does not contain natural rubber latex

## Ambu® Disposable Pressure Manometer

### Product Specification

Product item no.	322 003 000
Type	Ambu Disposable Pressure Manometer, 20 pcs per box sold individually
Measuring limits	0-60 cm H <sub>2</sub> O (hPa)
Measuring accuracy	+/- 2 cm H <sub>2</sub> O (hPa) tot 7,5, 10, 15, 20, 30 H <sub>2</sub> O (hPa) +/- 3 cm H <sub>2</sub> O (hPa) tot 40 H <sub>2</sub> O (hPa) +/- 5 cm H <sub>2</sub> O (hPa) tot 60 H <sub>2</sub> O (hPa)



## Notes

[illegible]

## Notes

[illegible]



# Ambu®

# Airway

Ambu Ltd.  
8 Burrell Road  
St Ives  
Cambridgeshire  
PE27 3LE

Tel 01480 498403  
Fax 01480 498405  
uksales@ambu.com  
www.ambu.co.uk

v3 06/2019

**Ambu**   
Ideas that work for life