

## New bronchoscopes, new possibilities

Say hello to the Ambu® aScope™ 4 Broncho. A new single-use, high quality bronchoscope and our latest proof of ideas that work for life. The aScope 4 Broncho is always available when you need it and sterile straight from the pack, offering improved patient safety and workflow. It comes in three sizes in one system at no additional cost. It is designed specifically for bronchoscopy and the ideal solution for a wide range of procedures.

"The ideal solution for a wide range of bronchoscopy procedures."

Instantly available, sterile straight from the pack for improved patient safety.

496 361 001 Vol. 2017/09

Experience new possibilities with aScope 4 Broncho.

Visit [ambu.co.uk](http://ambu.co.uk) today

**Ambu Ltd.** 8 Burrell Road St Ives Cambridgeshire PE27 3LE  
Tel 01480 498403 Fax 01480 498405 [uksales@ambu.com](mailto:uksales@ambu.com) [www.ambu.co.uk](http://www.ambu.co.uk)

Single-use. New possibilities.



# Patient safety you can rely on

### Improve your patients' safety

Reusable bronchoscopes are difficult to clean and disinfect due to their long narrow lumens and delicate materials. Even the most stringent cleaning procedures do not guarantee the complete removal of biofilm from a reusable scope. But the risk of device-related cross-contamination can be completely avoided with the sterile single-use aScope 4 Broncho. Simply use the scope and discard it.

The immediate access to the plug-and-play aScope 4 Broncho also removes the risks associated with waiting for an available broncho-scope to be delivered.

### Cost and time effective

The high quality and improved patient safety of the aScope 4 Broncho comes at a low investment cost. As a single-use system, you also avoid the reprocessing and repair costs associated with reusable bronchoscopes, because it is brand new, every time. The aScope 4 Broncho offers three readily available bronchoscopes in one system at no additional cost. So all your bronchoscopy needs are covered at a fixed budget per use.

The portable, ready-to-go system is easy to set up and intuitive to use, saving you valuable time that can be spent on your patients rather than on the equipment.



- aScope 4 Broncho Slim 3.8/1.2  
Bending: 180/180
- aScope 4 Broncho Regular 5.0/2.2  
Bending: 180/180
- aScope 4 Broncho Large 5.8/2.8  
Bending: 180/160

Three sizes in one system at no additional cost.


# High quality bronchoscopy

### High Quality view



aScope 4 Broncho connects to the Ambu® aView™, our portable, high-resolution monitor that is operational in seconds and intuitive to use.



aScope 4 Broncho is always there when you need it with no risk of cross contamination.



aScope 4 Broncho delivers clear visibility of near and distal bronchial segments.



Ergonomic design for comfortable grip and ease of use.



High bending angles and precise tip motion makes navigation easy.

### Clear and sharp images

The aScope 4 Broncho gives you clear visibility of the near and distal bronchial segments. The clear, sharp images and adaptive light control make identification of lung anatomy and navigation in the bronchial tree easy. Combined with a colour representation on the Ambu® aView monitor reflecting lung tissue colours, the aScope 4 Broncho offers high quality imaging for your bronchoscopy procedures.

### Smooth and easy navigation

The high bending angles and precise tip motion of aScope 4 Broncho allows you to easily navigate into all segments of the lung. The high quality handle is designed with a shaped, rippled surface to give you a comfortable grip throughout the procedure. And the lightweight and ergonomic design makes aScope 4 Broncho easy to use.

