



Ambu® aScope™ 4 Broncho

Guaranteed sterility
Ready when you are

Single-use – the smarter solution for enhanced patient safety

Guaranteed sterility

Even with the most stringent cleaning procedures, reusable endoscopes are very difficult to clean and complete sterility cannot be assured.

Single-use aScope 4 Broncho guarantees sterility straight from the pack, avoiding the risk of cross-contamination while improving safety, enabling you to focus on the best possible care for patients.





Spend time on your patients rather than the equipment

This portable, intuitive bronchoscopy solution is easy to set up, saving you valuable time and allowing you to concentrate on what really matters: your patients.

Ready when you are

The aScope 4 Broncho is ready when you are, reducing the risks and frustrations associated with waiting for an available endoscope.

The single-use bronchoscope requires zero handling or reprocessing after use, streamlining your workflow, optimising resources and boosting productivity.

Your hassle-free bronchoscopy solution

Intelligent, innovative design

The aScope 4 Broncho is more than a simple bronchoscope. It is a fully integrated solution, helping you to deliver the very best in patient care.

The ergonomic and lightweight design enhances the user experience throughout multiple procedures, while the scope's high bending angles and smooth and precise tip motion result in seamless access to all lung segments.

Single operator sampling solution

The aScope BronchoSampler is specially designed for use with the aScope 4 Broncho. This closed-loop system reduces the risk of sample loss and contamination, improving the safety of patients and carers.





Clear, crisp imagery

The aScope 4 Broncho connects to the Ambu® aView, our portable, high resolution monitor that delivers clear visibility of the bronchial segments.

The sharp imagery and adaptive light control allow for easy navigation of the bronchial tree, resulting in a quick and confident assessment of the bronchial tissue.

Three sizes for multiple procedures

To meet your needs across a wide range of bronchoscopy procedures, the aScope 4 Broncho comes in three sizes.

-  **aScope 4 Broncho Slim 3.8/1.2**
Bending: 180/180
-  **aScope 4 Broncho Regular 5.0/2.2**
Bending: 180/180
-  **aScope 4 Broncho Large 5.8/2.8**
Bending: 180/160



aScope BronchoSampler™

Breakthrough closed-loop sampling solution

Discover a world-first in bronchoscopy sampling with the breakthrough aScope BronchoSampler, a dedicated, compact sampling solution specifically designed for the aScope 4 Broncho.

A simplified sampling process

A one-of-a-kind, purpose-built system, the aScope BronchoSampler improves workflow efficiency by allowing a single clinician to easily perform bronchoscopic sampling as part of BAL or BW procedures – from system assembly to sample procurement.

aScope BronchoSampler also features a vacuum bypass that further simplifies procedures by enabling you to change between suction and sampling without a suction tube switch.

Fewer lost and contaminated samples

As a single-use, sterile, closed-loop system aScope BronchoSampler reduces the risk of sample loss and guarantees sample quality from start to finish. By minimising the risk of sample contamination or exposure, aScope BronchoSampler also supports accurate and timely diagnostics as well as the safety of both patients and clinicians.

Innovative, integrated system

aScope BronchoSampler is designed to integrate seamlessly with your aScope 4 Broncho. It comes with all the critical components that you need, allowing you to effectively and efficiently perform BAL and BW procedures. And because it's always ready when you are, aScope BronchoSampler minimises the significant risks to patients that can result from equipment and workflow delays.



Ambu Ltd.

8 Burrell Road
St Ives
Cambridgeshire
PE27 3LE
Tel 01480 498403
Fax 01480 498405
uksales@ambu.com

For more information

Visit ambu.co.uk

Item code	Product
aScope 4 Broncho flexible single-use bronchoscope	
476 001 000	aScope 4 Broncho Slim 3.8 / 1.2
477 001 000	aScope 4 Broncho Regular 5.0 / 2.2
478 001 000	aScope 4 Broncho Large 5.8 / 2.8
aView monitor	
405 002 000	aView monitor
aScope BronchoSampler	
480 200 000	aScope BronchoSampler, specimen sampling system
480 300 000	aScope BronchoSampler SC, sample container

About Ambu

Since 1937 Ambu® has been bringing breakthrough healthcare solutions to life within anaesthesia, patient monitoring and diagnostics, and emergency care. Today, we are dedicated to improving safety for millions of patients and healthcare professionals around the world through the advancement of single-use devices, including aScope – the world's first single-use flexible video scope. Our commitment to bringing new ideas, exceptional product performance and superior service to our customers has made Ambu® one of the most recognised medical companies in the world. Headquartered near Copenhagen, Denmark, Ambu® employs around 2,500 people across Europe, North America and Asia Pacific. For more information, please visit www.ambu.com